

Airbeds & Nebulizer

M-Care Air matters:

Operating Principles:

- Alternating pressure at sufficient rhythm avoids prolonged vascular compression which is likely to cause tissue hypoxia
- Alternating pressure adapts the matters to the patients and their position

General:

- Simplicity, reliability and easy to use
- Aluminium motor cover
- Quiet pump low noise, low vibration, low temperature, with knob easy adjust pressure
- Fixed hook easy fix on bed frame.
- For Patients who are likely to spend more than 15 hours per day in bed
- For Patients who are classed as medium to high risk of developing pressure ulcers.

Mattress Specification:

- Max. User Weight: 130Kgs.
- Mattress Dimensions: 200cm (L)* 90cm(W)*7.0cm(H)
- 130 Bubble Cells with end flaps
- PVC Thickness: 0.3mm
- Colour: Beige
- Mattress Weight: 2.2-2.5 kg

Compressor Characteristics:

- Power Supply: 220V/50Hz
- Air Output: 5-7 liter/min.
- Pressure Range: 50mmHg-110mmHg
- Cycle time: 12 min
- Size: 23.5(L)*12(W)*9.5(H)
- Weight: 1.4kg



General:

- Nebulizer is a drug delivery device used to administer medication in the form of a mist inhaled into the lungs.
- Nebulizers can be used for medical treatment of the upper and lower respiratory tracts.
- MINY-003 Compressor Nebulizer hordy, compact size and low noise features give a user friendly & easy respiratory treatment for you and child.

Specification:

- Input: 220V-50Hz
- Power Consumption: 180VA
- Medication Capacity: 6ml
- Average Nebulization Rate: Above 0.2ml/min
- Noise Level: Below 55dBA
- Compressor Pressure Range: 30-45Psi
- Operational Pressure Range: 12-14.5Psi (0.85-1.1 bar)
- Liter flow range: 8-10 lpm

M-Care Nebulizer



TECHNOLOGY OF KOREA



We support service...